

Installing Eye Tools

The *eye tools* program **cannot be run from the CD-Rom**. It must be installed on the hard drive of your computer to allow information and your own data to be entered for calculations.

To install, insert the *eye tools* CD in your CD drive..

- ◆ **If Autoplay is enabled**, the CD will start itself and begin the installation a few seconds after it is inserted in the CD drive.
- ◆ **If Autoplay is NOT enabled**, click on *START*, then *RUN*, then type D:setup.exe (where “D” is the letter of your CD drive).
- ◆ **After installation is complete**, start *eye tools* from the desktop icon. Then select *Activation* and enter your *User Name* and *Product Key*.

Document: 00102 Last edited: 20/5/05



Figure: The packaging for the software CD

Starting Eye Tools

To start the program, simply double-click on the *eye tools* icon on the desktop.




Figure: The desktop icon to launch eye tools