

Downloading New Data Files

Document: 00122 Last edited: 24/10/05

The *Load New Data* area is on the *Setups* screen (reached from the *Main Menu*). This allows a user to update the data that is used by the *eye tools* calculators.

✓ Tip:

- From time to time, new *Materials & Lenses* will become available and there will be additions or changes to the database of product names in *Brand Finder*.

To keep your program up-to-date, the latest data files are available for download from the *Eye Tools Support Knowledgebase*.

Downloading New Data Files

1. Go to the *Eye Tools Support Knowledgebase* at www.hanks.optom.com.au and check the dates of the most recent data files that are available:

Finder Update.USR
Calc'n Update.USR

2. Compare these to the dates of the files that you currently have installed in *eye tools*. (These are the *current data dates* displayed on the setups screen).
3. If more recent data files are available, click on "Download Now", then "Save".
4. Use the "Save In" window to navigate to your *eye tools* directory on your hard drive:

C:\Program Files\Eye Tools

5. Click on "Save", then answer "Yes" to replace the old file(s).

Loading the New Data Files

1. Start the *eye tools* program. Go to *Setups* and confirm that the date for the *New Data Available* has now changed. (If not, the files may not have been copied to the correct *eye tools* directory).
2. Click on the button *Load Now* to update your data to the more recent file(s). The button will then change to *Up to date*.
3. Click on the second *Load Now* button if both files are being updated.



Figure: The Eye Tools Support Knowledgebase at www.hanks.optom.com.au




Figure: The Load New Data area of the Setups screen shows the dates of the current data in the program